**size:60x72mm**

# User's Manual

## Infrared Ear Thermometer

Model **ET-101A**

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


Figure 12

## Introduction

The ET-101A thermometer is specifically designed for safe use in the eardrum. The thermometer measures infrared heat generated by the eardrum and surrounding tissues. It converts the measured heat into a temperature reading displayed on the LCD. When properly used, it will quickly assess your temperature in an accurate manner.

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
## Ear Thermometer Advantages

Infrared Ear Thermometer measures core body temperature, which is the temperature of a body's vital organs. (See Figure 1) Ear temperatures accurately reflect core body temperatures since the eardrum is in the cranial cavity and is not affected by the environment. The eardrum shares blood supply with the hypothalamus (part of the brain.) Therefore, body temperature changes are reflected sooner in the ear than they are in other parts of the body.

~Axillary (armpit) temperature measures skin temperature and therefore may not be a reliable indicator of core body temperature.

~Oral temperature can be influenced by drinking, eating, and breathing.

~Rectal temperature readings often lag behind changes in core body temperature. Across-contamination risk may exist as well.

## Body Temperature

The following table demonstrates how normal temperature ranges can vary by location of the body. Therefore measurements from different body locations should not be directly compared. Normal Ranges by Body Location:

Axillary	34.7 – 37.3°C	94.5 – 99.1°F
Oral	35.5 – 37.5°C	95.9 – 99.5°F
Ear	35.8 – 38.0°C	96.4 – 100.4°F

A person's normal temperature range tends to decrease with age.

Normal Ranges by Age:

0-2 Years	36.4 – 38.0 °C	97.5 – 100.4 °F
3-10 Years	36.1 – 37.8 °C	97.0 – 100.0 °F
11-65 Years	35.9 – 37.6 °C	96.6 – 99.7 °F
Over 65 Years	35.8 – 37.5 °C	96.4 – 99.5 °F

Normal Ranges vary from person to person and can fluctuate throughout the course of a day.

## Illustration For Use

1. Place a new probe cover on the sensor firmly; making sure the probe and the inner portion of the probe cover fit snugly. To ensure proper hygiene and to avoid contamination, a new probe cover should be applied for each user. If the probe cover is in good condition, it may be used several times by the same individual.
2. Press the On/Off Button or the Test Button. The LCD screen will appear as shown in Figure 2 (Note: The last recorded temperature reading will be present.) Wait for screen to change to setting shown in Figure 3. The thermometer is now ready for a new measurement. If there is an extraordinary, the Er1 or Er2 appears.
3. Place the probe snugly into the ear canal. Make sure the probe tip faces the eardrum as depicted in Figure 4. Press the Test Button. A long beep will sound approximately 2 seconds after the Test Button is activated. This signals the end of the process as results are displayed on the LCD screen.
4. The thermometer is ready for testing again once an ear mark displays on the LCD. A waiting period of 20 seconds between testing is recommended to ensure complete accuracy.
5. The thermometer will shut off automatically after 60 seconds of inactivity. To prolong battery life, press the On/Off Button to turn the unit off.

## Temperature Taking Hints

1. Temperatures of two eardrums may slightly differ. Always test the same ear drum when comparing results.
2. External factors may influence ear temperatures such as:
  - ˆ Lying on one ear or the other
  - ˆ Having one or both ears covered
  - ˆ Exposed to very hot or cold temperatures
  - ˆ Recent swimming or bathingUnder these circumstances, remain in the testing environment for 20 minutes or longer prior to temperature taking.
3. The ear must be free from obstructions such as excess earwax.

## Temperature Taking Hints

4. When using the thermometer on infants under age 1, pull the ear up making sure the sensor faces the eardrum. (See Figure 5)
5. When using the thermometer on individuals over the age of 1, pull the ear back making sure the sensor faces the eardrum. (See Figure 6)
6. Do not touch the probe tip. A cotton swab and rubbing alcohol may be used to clean the sensor tip.
7. If thermometer is stored in a significantly different environment than testing location, place it in the testing location for approximately 30 minutes prior to use.

## Memory Mode

1. The Memory Mode can be accessed two ways:
  - a) After the thermometer has been turned on and displays the screen depicted in Figure 3, press and hold the On/Off Button. The letter " M" will appear in the upper left corner of the display. (See Figure 7)
  - b) After testing, press and hold the On/Off Button for several seconds until the letter" M" appears in the upper left corner of the LCD.
2. The thermometer will automatically memorize the last 10 temperature readings. Each time the On/Off Button is pressed, the screen displays past readings that correspond with a number 0-9. The number 0 reflects the most recent reading while the number 9 reveals the oldest reading stored in memory. (See Figure 8)
3. After the LCD displays the last memorized reading(number 9) or press Test Button in the memory mode directly the thermometer will enter into testing mode.

## Selecting A Temperature Scale

1. Temperature readings are available in the Celsius (°C) or Fahrenheit (°F) scale.
2. With the unit off, press and hold the On/Off Button and continue holding until either C or F is present as shown in Figure 9 (this is the current setting.)
3. To change this setting to a different scale, release the On/Off Button and wait approximately 7 seconds until Figure 3 is present.
4. Press the On/Off Button once and the unit will turn off.
5. After the unit has shut off, press and hold the On/Off Button again until the new scale appears.
6. Release the On/Off Button and wait 7 seconds; the unit is now ready for testing in the newly chosen scale.

## Care And Cleaning

1. The probe tip is the most delicate part of the thermometer. To safely clean the tip, remove the probe cover first and gently wipe its surface with a cotton swab or soft cloth moistened with rubbing alcohol. Replace the probe cover only after the cleaned surface has completely dried.
2. Use a soft, dry cloth to clean the thermometer display and exterior.
3. The thermometer is not waterproof. Do not submerge the unit in water when cleaning.
4. Store thermometer and probe cover in a dry location free from dust and contamination. Do not place in direct sunlight.

## Battery Replacement

1. Battery needs replacing when the LCD displays the icon shown in Figure 10 and then turns off.
2. Slide battery cover down as shown in Figure 11.
3. Remove battery as shown in Figure 12.
4. Replace with one new battery (CR 2032.) Place the  $\sim + \sim$  side upward.

## Specifications

Accuracy:

$\pm 0.2^{\circ}\text{C} / \pm 0.4^{\circ}\text{F}$  ( $35.5^{\circ}\text{C} \sim 42.9^{\circ}\text{C}$  ;  $95.9^{\circ}\text{F} \sim 109.2^{\circ}\text{F}$ )

$\pm 0.3^{\circ}\text{C} / \pm 0.5^{\circ}\text{F}$  (under  $35.5^{\circ}\text{C} / 95.9^{\circ}\text{F}$ )

Response Time: Approximately 2 Seconds

Test Range:  $32.0^{\circ}\text{C} \sim 43.0^{\circ}\text{C}$  ( $89.6^{\circ}\text{F} \sim 109.4^{\circ}\text{F}$ )

Store and Transport Environment:

$-25^{\circ}\text{C} \sim 55^{\circ}\text{C}$  ( $-13^{\circ}\text{F} \sim 131^{\circ}\text{F}$ ),  $\leq 95\% \text{ RH}$ , non-condensing

Operating Environment:

$10^{\circ}\text{C} \sim 40^{\circ}\text{C}$  ( $50^{\circ}\text{F} \sim 104^{\circ}\text{F}$ ),  $\leq 95\% \text{ RH}$ , non-condensing

Dimension: 106x33x45mm

Weight: Approx. 34.5 grams

Battery: One 3.0V DC button battery (size CR2032) L

Battery Life: 2 years / 1000 measurements

Type B applied part

Direct Current


Correct Disposal of This Product

(Waste Electrical & Electronic Equipment)

## Troubleshooting

Situation	Solution	Error message
Button malfunctions	Remove old battery and replace with new one	No
Temperature results are not within typical human temperature range 32.0°C ~ 43.0 °C (89.6°F ~ 109. 4°F) Lo= too low;H = too high	Make sure the probe tip and lens filter are clean. Be sure thermometer is properly inserted into a clean ear canal. Contact seller if error message continues after several consecutive attempts.	
Ambient temperature is not within the allowed operating range 10°C ~ 40°C ; (50°F ~ 104°F)	Place thermometer in a room for 30 minutes with a temperature range between 10°C ~ 40°C (50°F ~ 104°F)	
Other error	Contact seller	

## Warranty

Thermometer complies with EN12470-5 Clinical Thermometer - Performance of Infrared Ear Thermometers, ASTM E-1965, IEC60601-1 and IEC60601-1-2. Manufacturer is ISO 9001 and ISO 13485 certified.

Thermometer is guaranteed for 1 Years from date of purchase. If thermometer does not function properly due to defective components or poor workmanship, we will repair or replace it free of charge. All components are covered by this warranty excluding the battery. This warranty does not cover damages due to improper use or handling.

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**Hangzhou Sejoy Electronics & Instruments Co.,Ltd**
**Shanghai International Holding Corp.
GmbH (Europe)**