

- Device classification: Class 1
- Duration: Temporary (rule 1)
- ANSM registration No: FR/CA01/0004344V

Applied rule (annex IX)

- Non-Invasive (regulation 5 # 1)
- Non active therapeutical device
- Non surgical
- Non active medical device
- Non active device for diagnosis

<i>General description</i>	Mobile infusion foot with hooks (2 or 4 plastic or stainless steel hooks) and 5-footed base, for carrying infusion bags or pouches. For the rods, simply see the height of the device without the feet.			
<i>Attachments</i>	Yes	Diagram + user instructions + parts list		
<i>Batch reference</i>	No			
<i>Connected device</i>	No			
<i>Single use device</i>	No			
<i>Sterile device</i>	No			
<i>Sterilisation</i>	No			
<i>Autoclavable</i>	No			
<i>Cleaning</i>	Yes	With a product that does not attack stainless steel or plastic in addition to soapy water or a disinfectant detergent.		
<i>Unsuitable substance</i>	Products that attack stainless steel, polyamide, and plastic			
<i>Component parts</i>	Yes	Rod: Supporting rod in AISI 304 quality stainless steel, polishing to a 220 grit level, dimensions 1000mm x Ø30mm x 1.2mm. Identical telescopic rod measuring 1200mm x Ø16mm x 1mm. Base: 5 feet in 20% fibreglass with polyamide caster wheels, 2 with clamps. Bottle holder: 2 or 4 nylon or stainless steel hooks.		
<i>Material</i>	Stainless steel rod, polyamide base, plastic or stainless steel hooks			
<i>Contains latex</i>	No			
<i>Labelling</i>	Yes	CE/MOD/NAME thermal printed at the top of the supporting rod		
<i>Packaging</i>	Yes	Cardboard with plastic sleeve for the rods		
<i>Special packaging</i>	No			
<i>Marking on medical device</i>	No			
<i>References</i>	45 400 05	45 404 05	45 410 05	45 410 06
<i>Designation</i>	2 plastic hooks	4 plastic hooks	2 stainless steel hooks	4 stainless steel hooks
<i>Rod references</i>	4540005T	4540405T	4541005T	4541006T
<i>Dimensions (cm)</i>	Appendix 1			
<i>Net weight</i>	2900	2910	2950	3000
<i>Temperature resistance (C°)</i>	Max = + 40, Min = -5			
<i>Load weight (kg)</i>	Maximum: 2 per hook			
<i>Usage</i>	For hanging infusion bags or pouches.			
<i>Instruction for use</i>	Attach the infusion bottles or pouches to the hooks and adjust to the desired height using the telescopic rod, then tighten using the knob. If necessary, lock the stand in place by applying the wheel clamps (2).			
<i>Precautions for use</i>	The hooks can hold a maximum weight of 2kg. The maximum height is 2100mm. Ensure the floor surface is flat or at a gradient of less than 5.7° (10%) in order to prevent toppling. Never lift with pouches attached. Always extend the base by using the rod, holding it 25mm from the locking screw.			
<i>Recycling/Recovery</i>	No/In metal waste containers.			
<i>WEEE/REACH</i>	No/No			
<i>Safety data</i>	R-Phrases = None Combination phrases = None			
<i>Applicable standards</i>	NF EN ISO 9001 : 2015 ; NF EN ISO 13485 : 2012 NF EN 1041 : 2008 + A1 : 2013 NF EN ISO 14971 : 2013 NF EN ISO 15223-1 : 2012 NF ISO 15223-2 : 2010			

ANNEX 1:

To assemble, insert the rod with the black tip into the 5-footed base and fix in place using the bolt, placing the washer underneath the frame. Clean the two rods with a cleaning product before using for the first time. 2 lockable wheels are provided.

REF 45 400 05 = 2 plastic hooks (1)

REF 45 410 05 = 2 stainless steel hooks (3)

REF 45 404 05 = 4 plastic hooks (2)

REF 45 410 06 = 4 stainless steel hooks (4)

