

<i>General description</i>	Dermatograph marker	
<i>Enclosed pieces</i>	Yes	Picture.
<i>Lot reference</i>	No	
<i>Connected device</i>	No	
<i>Single Use Device</i>	No	
<i>Dispositif stérile</i>	No	
<i>Sterilization</i>	No	
<i>Autoclaveable</i>	No	
<i>Cleaning</i>	No	Can not be cleaned
<i>Unfitting substance</i>	Product hurtful to polyethylene and polypropylene	
<i>Constituent pieces</i>	No	
<i>Content</i>	Cap made of polyethylene, Body in polypropylene Washable ink without solvent mainly made of water, cosmetic grade colour dye and food grade preservative.	
<i>Labelling</i>	No	
<i>Packing</i>	No	
<i>Special storage</i>	Yes	Dry and clean
<i>Marking on the Device</i>	Yes	COMED
<i>References</i>	29 160 40	29 160 45
<i>Designation</i>	Red Marker	Black marker
<i>Dimensions Length x Ø (mm)</i>	138 x 12.5	
<i>Use life</i>	1 year after first opening	
<i>Temperature resistance (C°)</i>	-5 to 35°C	
<i>Netto weight (g)</i>	8	
<i>Colours</i>	Red	Black
<i>Use</i>	Device used to mark on skin, marker.	
<i>Instructions for use</i>	Take a marker and mark on desired space. To wash patient skin, use water and soap. Repeat process as many times as needed. Number of time will depend on patient skin type. The marker does not fit on betadine.	
<i>Use care</i>	Correctly close cap after each use.	
<i>Recycling/Recovery</i>	In medical waste containers.	
<i>PHTALATES</i>	No/ NC	
<i>DEEE</i>	No/ No	
<i>Security data</i>	Phrases R = Blank      Phrases combination = Blank	