

Cotton balls



COMPOSITION	<ul style="list-style-type: none"> • 100% bleached absorbent cotton or cotton/viscose mix • H2O2 bleached (non chlorine bleached, free of any optical brightener) • Meeting the European Pharmacopoeia Standards & Oeko-tex 100 certified • In accordance with REACH European standards • Bioburden analysis on request
COTTON BALLS	<ul style="list-style-type: none"> • Different sizes available: from 0,3 to 5 g • Used for absorption of exudates, peripheral cleaning wounds and supporting the application of products on the skin + laboratory plugs • EO (Ethylene Oxide) sterilizable • Shelf life : 5 years from production date

Standard sizes and packaging

Reference	Description	Balls/Bag	Weight/Bag	Bags/box	Boxes/Pallet
1500B03	1500 balls 0.3 g	1500	(450g)	12	16
100B06COT*	100 balls 0.6g	100	(60g)	25	36
100B06COTC*	100 coloured balls 0.6g	100	(60g)	25	36
500B06	500 balls 0.6-0.7g	500	(325g)	16	16
700B06	700 balls 0.6- 0.7g	700	(455g)	16	16
500G03	500g of balls 0.3g	(1666)	500g	10	16
500G06	500g of balls 0.6g	(833)	500g	10	16
500G09	500g of balls 0.9g	(555)	500g	10	16
5KG06	5 kg of balls 0.6-0.7g	(7692)	5000g	1	16
5KG2	5 kg of balls 2g	(2500)	5000g	1	16
5KG3	5 kg of balls 3g	(1666)	5000g	1	16
5KG5	5 kg of balls 5g	(1000)	5000g	1	16

Comments :

Class I medical device, single use product
 Packed in polyethylene bags to protect from dust and humidity
 The goods can also be delivered in bulk
 Other packing and packaging available on request
 * Available in Cotoneige brand
 Export box sizes: 60x40x60 cm (16) or 59x39x24 cm (36)
 General tolerance: +/- 10%
 Should be stored in a clean and dry place, out of direct sunlight

We are ISO 13485 certified and our products are CE marked



ALAN & CO

Avenue André Ernst 17 – 4800 Verviers – Belgium
 Tel : +32 87 31 10 59 – Fax : +32 87 31 61 69
 www.alan.be – info@alan.be