

COMED

STRASBOURG

**PULSE OXIMETER POUCH**  
 Ref. No.: 23 200 47 / 23 200 52  
 TCE2

TECHNICAL AND  
SAFETY SPECIFICATIONS

forming part of the technical file

<i>General description</i>	Pulse oximeter pouches	
<i>Attachments</i>	Yes	Photo
<i>Serial reference</i>	Yes	
<i>Connected device</i>	No	
<i>Single use device</i>	No	
<i>Sterile device</i>	No	
<i>Sterilisation</i>	No	
<i>Autoclavable</i>	No	
<i>Cleaning</i>	Yes	Machine wash
<i>Unsuitable substance</i>	Products that attack synthetics and cotton	
<i>Component parts</i>	No	
<i>Material</i>	Synthetic fabric and cotton	
<i>Contains latex</i>	No	
<i>Tube</i>		
<i>Standard manufacturing</i>	N/A	
<i>Labelling</i>	No	
<i>Packaging</i>	Yes	In a bag
<i>Special packaging</i>	No	
<i>Marking on medical device</i>	No	
<i>References</i>	23 200 47	23 200 52
<i>Designation</i>	Green pouch	Orange pouch
<i>Dimensions L x W x H (mm)</i>	90 x 63 x 35	
<i>Temperature resistance (C°)</i>	- 40° to 95°	
<i>Net weight (g)</i>	20	
<i>Colours</i>	green	orange
<i>Usage</i>	For carrying and protecting a pulse oximeter	
<i>Instruction for use</i>	Undo Velcro, insert pulse oximeter and reattach Velcro. May be attached to a belt or hook.	
<i>Precautions for use</i>	No	
<i>Recycling/ Recovery</i>	Yes, in household waste	
<i>WEEE/REACH</i>	No/NC	
<i>PVC/PHTHALATES</i>	No/No	
<i>Safety data</i>	R-Phrases = None    Combination phrases = None	