

Zig-Zag cotton wool



COMPOSITION	<ul style="list-style-type: none"> • 100 % bleached absorbent cotton or cotton/viscose mix • H2O2 bleached (non chlorine bleached, free of any optical brightener) • Complying with the European Pharmacopoeia Standards & Oeko-tex 100 certified • In accordance with REACH European standards • Bioburden analysis on request
ZIG-ZAG COTTON WOOL	<ul style="list-style-type: none"> • Used for absorption of exudates, cleaning peripheral wounds and supporting the application of products on the skin • EO (Ethylene Oxide) sterilizable • Shelf life: 5 years from production date

Packaging

Reference	Weight	Bags/Box	Boxes/Pallet
OU50G*	50 g	150	16
OU100G*	100 g	80	16
OU200G*	200 g	50	16
OU250G*	250 g	24	30
OU500G*	500 g	25	16

Comments :

Class I medical device, single use product
 Packed in polyethylene bags (with strings) to protect from dust and humidity
 Pre-cut (except for OU500G) on request
 Other packing or packaging available on request
 *Available in Cotoneige brand
 ** Minimum order quantity of 15 pallets or according to the stock
 Export box sizes: 60x40x60 cm (per 16) and 40x40x45 cm (per 30)
 General tolerance: +/- 10%
 Should be stored in a dry and clean place, out of direct sunlight

We are ISO 13485 certified and our products are CE marked

