

Medical range

CarePulp® Receptacles



Kidney dish

- Light and strong, the kidney dish is a universal receptacle for receiving soiled dressings and medical waste.
- The shape of the dish allows it to be held against the patient's body.

Technical data

Capacity	Package dimensions in mm (L x W x H)	Qty/ carton	Package weight in kg	Colours available
750ml	590x400x265	320	7	grey

Small pallet		
Height mm	Number of cartons per pallet	Qty/ pallet
1 475	4 x 5	6 400

Large pallet		
Height mm	Number of cartons per pallet	Qty/ pallet
2 800	4 x 10	12 800

Average weight per piece: 22 g

Other packaging and pallet configurations on request

Pallet dimensions (L x l): 800 x 1 200 mm

Disposable pulp products for use in healthcare

Conforms to British Standard PAS 29 : 1999

All products are waterproof for a minimum of 4 hours

Shelf life: 3 years

We recommend that receptacles are emptied on a timely basis to avoid the possibility of leakage after this period

The product is totally maceratable. Contact us for any information