

<i>General description</i>	Plasticized metal hook or not. Handheld system rounded or square pockets for pre-hole urine or other with loop for passing.	
<i>Enclosed documentation</i>	Yes	Photo
<i>Batch reference number</i>	No	
<i>Connected device</i>	No	
<i>Single Use Device</i>	No	
<i>Sterile Device</i>	No	
<i>Stérilization</i>	Yes	For the references 45 500 50/ 45 500 51, autoclave
<i>Autoclavable</i>	Yes	134°C
<i>Cleaning</i>	Yes	Soapy water or detergent which does not attack chrome, steel and plastic
<i>Unsuitable substances</i>	products Attacker the steel wire or plastic	
<i>Product components</i>	No	
<i>Materials</i>	Made of wire steel (45 500 50/51) or plastic (45 500 50R/51R)	
<i>Labeling</i>	No	
<i>Packing</i>	Yes	In an individual plastic bag
<i>Particular storage</i>	Yes	dry, clean, temperature between -40 and 95 ° C
<i>Marking on the Medical Device</i>	No	
<i>Références + extension</i>	45 500 51/51R	45 500 50/50R
<i>Désignation</i>	rounded hook	square hook
<i>Dimensions Length x Width x depthx diameter (mm)</i>	153x110x61x3.05	157x95x52x4.05
<i>Capacity (kg)</i>	2	2
<i>Temperature resistance (C°)</i>	-40°C at + 95°C	
<i>Net Weight (g)+- 1%</i>	44	60
<i>Colors</i>	Chrome or white	
<i>Use</i>	Apparatus used to support the pockets especially the urine pockets at the foot of the care beds.	
<i>Use care</i>	Well hang the support before determining if support flange for the hook is round or square. Choosing the appropriate support Wash after each use. if the plastic coating is damaged, risk of rust in this case no longer use.	
<i>Recycling</i>	Yes	
<i>Recovery</i>	In a medical waste bin contaminated	
<i>DEEE/REACH</i>	No/No	
<i>Phthalates/PVC</i>	No/No	
<i>Security data</i>	R Phrases = None	Phrases combination=None