

# Medical range

## CarePulp® Urinal



Male round urinal

- Traditional style, easy for staff to transport.
- Enables hygienic, independent toileting.

### Technical data

Capacity	Package dimensions in mm (L x W x H)	Qty/ carton	Package weight in kg	Colours available
875ml	945x255x885	100	5	grey

Small pallet		
Height mm	Number of cartons per pallet	Qty/ pallet
1 035	5	500

Large pallet		
Height mm	Number of cartons per pallet	Qty/ pallet
2 805	5 x 3	1 500

Average weight per piece: 50 g

Other packaging and pallet configurations on request

Pallet dimensions (L x l): 1 200 x 1 000 mm

Disposable pulp products for use in healthcare

Conforms to British Standard PAS 29 : 1999

All products are waterproof for a minimum of 4 hours

Shelf life: 3 years

We recommend that receptacles are emptied on a timely basis to avoid the possibility of leakage after this period

The product is totally maceratable. Contact us for any information