



<i>General description</i>	Pill crusher	
<i>Attachments</i>	Yes	Photo
<i>Batch reference</i>	Yes	
<i>ACL Code</i>	No	
<i>Connected device</i>	No	
<i>Single use device</i>	No	
<i>Sterile device</i>	No	
<i>Sterilisation</i>	Yes	In ethylene oxide
<i>Autoclavable</i>	No	
<i>Cleaning</i>	Yes	Water + household cleaner
<i>Unsuitable substances</i>	Substances that attack plastic (ABS, PS)	
<i>Component parts</i>	No	
<i>Material</i>	ABS, PS	
<i>Labelling</i>	Yes	Reference, batch, COMED, denomination in 7 languages
<i>Packaging</i>	Yes	Individual
<i>Special packaging</i>	Yes	In a clean, dry place away from sunlight
<i>Duration of use</i>	Short	
<i>Marking on medical device</i>	No	
<i>References + extensions</i>	16 200 52	
<i>Designation</i>	Ergonomic pill crusher	
<i>Dimensions L x W x H (mm) +- 3mm</i>	68 x 50 x 48	
<i>Net weight (g)+1%</i>	49	
<i>Colours</i>	Blue and transparent	
<i>Temperature resistance (°C)</i>	-10 to +60°C	
<i>Usage</i>	Device for crushing pills to facilitate ingestion.	
<i>Instruction for use</i>	Place the pill in the crusher and slowly screw the top part to crush. To make crushing easier and avoid damaging the crusher, it is highly recommended that the pills be placed sideways on and not flat as this will bring the threaded surface closer and make crushing easier.	
<i>Precautions for use and/or single use warning</i>	Thoroughly clean the crusher before and after use to prevent two different medications from mixing.	
<i>Recycling</i>	No	
<i>Recovery</i>	In household waste containers	
<i>Latex/Antistatic</i>	No/No	
<i>WEEE/REACH</i>	No/No	
<i>Phthalates/PVC</i>	No/No	



Technical and Safety Specification

Pill crusher

REF. 16 200 52

TCE3 - Non-medical device

Created: 26.11.2012

Update: 21.02.2018

Revision letter: c

Page 2 / 2

RoHS	NC
Safety data	R-Phrases = None Combination phrases = None

Editor: Marion KLEIN
Role: QA
Initials: MK

COMED SAS
8 rue Louise Michel
67200 STRASBOURG - France
www.comed.fr