



|  |  |  |
|--|--|--|
| <i>General description</i>                           | Weekly pill organiser  |  |
| <i>Attachments</i>                                   | Yes  | Photo  |
| <i>Batch reference</i>                               | Yes  |  |
| <i>ACL Code</i>                                      | No   |  |
| <i>Connected device</i>                              | No   |  |
| <i>Single use device</i>                             | No   |  |
| <i>Sterile device</i>                                | No   |  |
| <i>Sterilisation</i>                                 | No   |  |
| <i>Autoclavable</i>                                  | No   |  |
| <i>Cleaning</i>                                      | Yes  | Water + household cleaner  |
| <i>Unsuitable substances</i>                         | Substances that attack plastic (ABS, PS, PP)   |  |
| <i>Component parts</i>                               | No   |  |
| <i>Material</i>                                      | ABS, PS, and PP  |  |
| <i>Labelling</i>                                     | Yes  | On the box: Reference, batch, COMED, denomination in 7 languages |
| <i>Packaging</i>                                     | Yes  | Individual   |
| <i>Special packaging</i>                             | Yes  | In a clean, dry place away from sunlight                         |
| <i>Duration of use</i>                               | Temporary  |  |
| <i>Marking on medical device</i>                     | Yes  | Comed logo + days of the week + "Morning/Midday/Evening/Night"   |
| <i>References + extensions</i>                       | 16 200 56  |  |
| <i>Designation</i>                                   | Weekly pill organiser  |  |
| <i>Dimensions L x W x H (mm) +- 3mm</i>              | 115 x 43 x 128   |  |
| <i>Net weight (g)+1%</i>                             | 171  |  |
| <i>Colours</i>                                       | White and blue   |  |
| <i>Temperature resistance (°C)</i>                   | -10 to +40°C   |  |
| <i>Usage</i>   | Item for preparing and organising pills, up to 4 doses per day for seven days.   |  |
| <i>Instruction for use</i>                           | Remove the pill boxes from their holder via the opening at the bottom of the holder and replace them by inserting them into the opening at the top of the holder. Slide open the pill box lid, place the pills in the corresponding slots, then slide the lid back on to close the pill box. |  |
| <i>Precautions for use and/or single use warning</i> | The boxes should only be loaded by a knowledgeable person as there is a risk of administering the wrong treatment if not used correctly.   |  |
| <i>Recycling</i>                                     | No   |  |
| <i>Recovery</i>                                      | In household waste containers  |  |



Technical and Safety Specification  
Weekly pill organiser  
REF. 16 200 56  
TCE3 - Non-medical device

Created: 18.07.2017  
Update: 18.07.2017  
Revision letter: a

Page 2 / 2

|                                     |  |
|-------------------------------------|--|
| <i>Latex/Antistatic/Bisphenol A</i> | No/No/No                                       |
| <i>WEEE/REACH</i>                   | No/No  |
| <i>Phthalates/PVC</i>               | No/No  |
| <i>RoHS</i>                         | NC   |
| <i>Safety data</i>                  | R-Phrases = None    Combination phrases = None |

Editor: Marion KLEIN  
Role: QA  
Initials: MK

COMED SAS  
8 rue Louise Michel  
67200 STRASBOURG - France  
[www.comed.fr](http://www.comed.fr)