



<i>General description</i>	Toothbrush impregnated with toothpaste	
<i>Attachments</i>	Yes	Photo
<i>Batch reference</i>	Yes	
<i>ACL Code</i>	No	
<i>Connected device</i>	No	
<i>Single use device</i>	No	
<i>Sterile device</i>	No	
<i>Sterilisation</i>	No	
<i>Autoclavable</i>	No	
<i>Cleaning</i>	No	
<i>Unsuitable substances</i>	Substances that attack plastic (PS, PP)	
<i>Component parts</i>	No	
<i>Material</i>	PS (handle) and PP (brush) Toothpaste ingredients: water, sorbitol, glycerol, K12, DL-menthol, flavouring, thickener, softener, carbomer, butanamide	
<i>Labelling</i>	Yes	On the box: Reference, batch, COMED, denomination in 7 languages
<i>Packaging</i>	Yes	Individual then 100 units per box
<i>Special packaging</i>	Yes	In a clean, dry place away from sunlight
<i>Duration of use</i>	Short	
<i>Marking on medical device</i>	No	
<i>References + extensions</i>	16 700 00	
<i>Designation</i>	Toothbrush impregnated with toothpaste	
<i>Dimensions L x W x H (mm) +- 3mm</i>	155 x 12 x 14	
<i>Net weight (g) +- 1%</i>	5	
<i>Colours</i>	White and blue	
<i>Temperature resistance (°C)</i>	-10 to +40°C	
<i>Usage</i>	Oral hygiene device for cleaning teeth	
<i>Instruction for use</i>	Remove the toothbrush from its plastic sheath, moisten the brush with a little water if necessary, then proceed to brush teeth.	
<i>Precautions for use and/or single use warning</i>	Single-use item, discard after use. Risk of contamination from shared use.	
<i>Recycling</i>	No	
<i>Recovery</i>	In household waste containers	



Technical and Safety Specification  
Toothbrush impregnated with toothpaste  
REF. 16 700 00  
TCE3 - Non-medical device

Created: 19.07.2017  
Update: 19.07.2017  
Revision letter: a

Page 2/2

<i>Latex/Antistatic/Bisphenol A</i>	No/No/No
<i>WEEE/REACH</i>	No/No
<i>Phthalates/PVC</i>	No/No
<i>RoHS</i>	NC
<i>Safety data</i>	R-Phrases = None    Combination phrases = None

Editor: Marion KLEIN  
Role: QA  
Initials: MK

COMED SAS  
8 rue Louise Michel  
67200 STRASBOURG - France  
[www.comed.fr](http://www.comed.fr)